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**Franck-Hertz Supply Unit**  
**Socket for FH tube, Hg with DIN Plug**  
**Ne-FH Connector Cable, 7-pin**

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The devices described in the following permit investigation of the discrete energy absorption of atoms for collisions with free electrons (electron collision experiment of Franck and Hertz). Experiments can be carried out both with mercury atoms in a heated Franck-Hertz tube, Hg (555 85 in electric oven 555 81 for 230 V or 555 82 for 115 V) and with neon atoms (neon Franck-Hertz tube 555 870 with socket 555 871).

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

\*) ESD = electrostatic discharge

## 1 Safety notes

- Do not switch on the supply unit until after connecting the Franck- Hertz tube.
- When using the Hg tube:  
Ground the copper cylinder to the yellow and green safety socket via the lead ⓑ (Fig. 3)!  
Connect the electric oven before switching on the supply unit, as mains voltage can be present at the output sockets Ⓢ on the rear panel of the apparatus.  
Be sure to insert the temperature sensor all the way into the corresponding blind hole ⓐ (Fig. 2) of the copper tube to prevent overheating of the Hg tube.  
A blinking display means that there is an error in the temperature measurement setup (see 3.3.1)  
If the collector current in the Hg tube suddenly rises, immediately switch the acceleration voltage  $U_2$  to 0 V. The gas discharge can be seen in the tube as a blue luminance. (s. 3.3.4)
- When using the Ne tube:  
Disturbances in the curve due to mains interference, particularly in experiments using the oscilloscope, can be minimized by further screening and by grounding conductive parts in the vicinity.
- The sensitive electronics of this apparatus can be impaired or damaged by electrostatic discharge (ESD).  
Measures:  
– choose the experiment area so that no electrostatic charges can build up on either the experimenter or the devices, or  
– ensure that such charges are discharged safely, e.g. via a grounded metal rod (connection rod, Cat. No. 532 16).

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Fig. 1.2  
Belegung der DIN-Buchse (1.1)  
Pin assignment of DIN socket (1.1)

#### ① Connection field

(1.1) Female DIN connector for connecting the tube socket of the Franck-Hertz tube, Hg with DIN plug (555 861) or for the neon FH connecting cable, 7-pin (555 872), with automatic recognition of the Ne tube  
Pin assignment see Fig. 1.2

(1.2) 4-mm sockets for connection socket of Franck-Hertz tube (Hg) 555 86

Designations: f, f<sub>k</sub>: Cathode heating element  
f<sub>k</sub>: Cathode  
g<sub>1</sub>: Emission grid  
g<sub>2</sub>: Anode grid  
A: Collector  
S: Discharge pins

(1.3) Connector slot for connection to CASSY interface (e.g. from 524 007)

#### ② Parameter and display field

- (2.1) DIN socket for connecting the temperature sensor  
(2.2) Screwdriver potentiometer for setting the temperature of the electric oven for the Franck-Hertz tube, Hg:  
Range: approx. 140°C to 210°C  
Preset value: 180°C  
(2.3) Potentiometer for setting the negative field voltage  $U_3$ : 0 ... 10 V  
(2.4) Potentiometer for setting the driving potential or control voltage  $U_1$ : 0...5V  
(2.5) Selector switch for parameters and measured values:  
Set temperature:  $\vartheta_s$   
Actual temperature (measured value):  $\vartheta$   
Driving potential:  $U_1$   
Negative field voltage:  $U_3$   
Acceleration voltage:  $U_2$   
Measured value:  $I$  (collector current)  
(2.6) Digital display of the parameter or measured value set with (2.5):  
Display in °C, V or nA; negative field voltage  $U_3$  is a negative value

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(2.7) Analog output (4-mm sockets) for voltmeter (pointer meter, zero point left) for displaying the quantity selected with (2.5) resp. the voltage proportional to that quantity:

The following applies for:  $\vartheta_s$  and  $\vartheta$ :  $100^\circ\text{C} \cong 1\text{ V}$   
 $U_1$ :  $1\text{ V} \cong 1\text{ V}$   
 $U_3$ :  $-1\text{ V} \cong 1\text{ V}$   
 $U_2$ :  $10\text{ V} \cong 1\text{ V}$   
 $I$ :  $1\text{ nA} \cong 1\text{ V}$

③ Control field

(3.1) Mode switch:

/// fast sawtooth signal for oscilloscope  
 ↘ slow ramp signal for recorder  
 RESET resets acceleration voltage to 0 V  
 MAN. manual setting with (3.2) for point-by-point recording  
 CASSY measurement with CASSY via (1.3)

(3.2) Potentiometer for manually setting the acceleration voltage  $U_2$ :

Hg: 0...30 V  
 Ne: 0...80 V

(3.3) Schematic circuit diagram of the experiment setup

(3.4) LEDs for indicating the type of tube connected:

Hg: green briefly, immediately after switching on the supply unit;  
 red when the oven is heating up;  
 green when the system is in the control range (set oven temperature)  
 Ne: green, automatic recognition via DIN socket (1.2)

(3.5) 4-mm output sockets for acceleration voltage  $U_2 = 10 \cdot U$

④ Collector current output field with 4-mm output sockets for the voltage  $U_A$  proportional to the collector current  $I_A$ :  $1\text{ V} \cong 1\text{ nA}$

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⑤ Mains connection with primary-side fuse  
 Connection voltage: 230 V (or 115 V, according to label), 50/60 Hz  
 Fuse: T 2.5 at 230 V and 115 V (698 19)  
 Power consumption: max. 30 W  
 approx. 230 W with electric oven

⑥ Mains switch

⑦ Screwdriver potentiometer for setting the cathode heating:  
 Voltage range: approx. 5.9 V... 6.7 V, default setting: 6.3 V

⑧ 4-mm safety sockets for connecting the electric oven (555 81 or 82) for the Franck-Hertz tube, Hg

Dimensions: 21 cm \* 30 cm \* 23 cm

Weight: 3 kg

Included in scope of supply (not shown): 14-conductor connector cable for connecting a CASSY interface to the supply unit via the connector slot (1.3).

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## 2.2 Socket for Franck-Hertz tube, Hg with DIN plug (555 861)

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Fig. 2 Socket with copper tube

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For connecting the Franck-Hertz tube, Hg (555 85) to the DIN socket (1.2) of the supply unit.

Ceramic socket with 6-conductor cable; conductor for collector current and DIN plug screened. With built-in series resistor for cathode heater.

Copper tube for matching the Franck-Hertz tube, Hg, to the diameter of the electric oven (555 81) and for screening the Hg tube from interference fields by grounding via a copper lead **b** with 4-mm plug; with blind hole **a** for the measuring tip of a temperature sensor.

### 2.3 Neon FH connector cable, 7-pin (555 872)

For connecting the socket for neon Franck-Hertz tube (555 871) with attached neon Franck-Hertz tube (555 870) via the DIN socket (1.2) on the supply unit.

6-conductor cable; conductor for collector current and DIN plug screened.

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### 3 Operation

#### 3.1 Equipment required

##### 3.1.1 Experiment procedure with the Franck-Hertz tube, Hg (555 85):

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
Fig. 3  
Connection of the Franck-Hertz tube, Hg with oven and temperature sensor. Insert temperature sensor in the blind hole (a) (see also Fig. 2) all the way; ground the copper tube to the yellow-green safety socket using the copper lead (b).

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1 Franck-Hertz tube, Hg	555 85
1 Connection socket for FH tube (Hg) with DIN plug	555 861
oder	
1 Connection socket for FH tube (Hg) with 4-mm plugs (see inset diagram Fig. 3)	555 86
1 Electric oven for 230 V AC	555 81
or for 115 V AC	555 82
1 Temperature sensor	666 193

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**3.1.2 Experiment procedure with the neon Franck-Hertz tube (555 870):**

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- Connecting the Ne-FH tube
- 1 Ne-FH connector cable 555 872
  - 1 Holder for Ne-FH tube 555 871
  - 1 Neon Franck-Hertz tube 555 870

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**3.2 Measuring equipment**

**3.2.1 Manual point-wise recording of a curve:**

You do not absolutely require an additional measuring device to record a curve point by point. Voltage and current values can be read from the display (2.6) of the supply unit by switching the selector switch (2.5) between  $U_2$  and  $I$ .

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Point-by-point recording of a curve with a pointer meter, measuring range 3/10 V DC (e.g. 531 94 or 531 911) Mode switch (3.1) is set to MAN. The acceleration voltage is set using the potentiometer (3.2). Measured values are recorded by switching the selector switch (2.5) between positions  $U_2$  and  $I$

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Point-by-point recording of a curve with two pointer meters, measuring range 3/10 V DC (e.g. 531 94 or 531 911) Mode switch (3.1) set to MAN. or  $\nearrow$  Set the voltage  $U_2$  with (3.2)

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### 3.2.2 Automatic recording with an XY-recorder:


Fig. 6

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Recording the curve with an XY-recorder (e.g. 575 622) Set the mode switch (3.1) to  $\nearrow$  so that the acceleration voltage increases continuously. Raise the pen of the XY-recorder before switching to RESET.

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Fig. 7

Oscilloscope display of the curve Set mode switch (3.1) to  $\swarrow$  so that the acceleration voltage is applied to the FH tube as a sawtooth signal.

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### 3.2.4 Computer-supported measured-value recording with CASSY:


Fig. 8

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Fig. 8

Computer-supported measured value recording with CASSY In computer-supported recording of the Franck-Hertz curve, the acceleration voltage is controlled by the software via the CASSY interface (set mode switch (3.1) to CASSY)

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- 1 CASSYpackE 524 007
- 1 Software "Universal Data Acquisition" 525 032
- You can use the measuring example supplied on the disk to record measured values for your own experiments

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### 3.3 Experiment procedure

Set up the experiment as shown in Figs. 3 - 8 and switch on the supply unit.

#### 3.3.1 Setting the temperature (see Fig. 1.1, Fig. 3)

Only required for experiments with the Franck-Hertz tube, Hg. The neon Franck-Hertz tube operates at room temperature, and is ready for operation after approx. 1 min (cathode warm-up time).

Connect the electric oven to ⑧; insert the temperature sensor in the blind hole ⑨ and connect it to the DIN socket (2.1). Then switch on the supply unit with ⑥: after a few seconds, the LED (3.4) for the Hg tube lights up red (oven heating up).

Wait until the electric oven reaches the operating temperature of the Hg tube (approx. 10 to 15 min.): The LED indicator (3.4) for the Hg tube turns green.

**Attention!** A blinking display (2.6) means that there is an error in the temperature measuring setup:

- The temperature sensor is not connected to the DIN socket (2.1). The display starts to blink immediately when the supply unit is switched on.
- The temperature sensor is not inserted in the blind hole of the copper cylinder. Blinking starts some minutes after the supply unit is switched on: the oven may have heated up somewhat!
- The temperature sensor has fallen out of the blind hole of the copper cylinder after the operating temperature has been reached. Blinking starts a few seconds after the sensor falls out.

#### 3.3.2 Setting the measuring instruments

- Acceleration voltage  $U_2$ :  
10 V  $\pm$  1V at output (3.5) or (2.7)

	Hg	Ne
Voltmeter:	0...3 V DC	0...10 V DC
XY-recorder, X-axis:	Var. approx. 0,1 V/cm	Var. approx. 0,3 V/cm
Oscilloscope, X-axis:	Cal. 0,5 V/cm	Cal. 1 V/cm
Collector current $I$ :	1 nA $\pm$ 1 V at output ④	
Voltmeter:	0...15 V	
XY-recorder, Y-axis:	Var. approx. 1 V/cm	
Oscilloscope, Y-Axis:	Cal. approx. 2 V/cm	

#### 3.3.3 Setting the tube parameters

- Experiments without oscilloscope (Fig. 5.1/ 2, Fig. 6): Set mode switch (3.1) to MAN.: you can vary the acceleration voltage  $U_2$  using the potentiometer (3.2).

Setting the parameter voltages  $U_1$  and  $U_3$ :

Typical values: Hg tube:  $U_1 = 1.5$  V  $U_3 = 1.5$  V  
Ne tube:  $U_1 = 3$  V  $U_3 = 8$  V

Slowly increase the acceleration voltage and observe the curve of the collector current  $I$  (resp. the voltage  $U_A$  proportional to it).

- Experiments with oscilloscope (Fig. 7): Set mode switch (3.1) to  $\infty$ : The acceleration voltage  $U_2$  is automatically run up to 30 V resp. 80 V in a sawtooth signal. On the oscilloscope screen, you can observe a curve like the one shown in Fig. 9.1.

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Display of curve on screen without parameter voltages  $U_1$  and  $U_3$   
The shape of the curve is influenced by discharge currents caused by the capacitances of the cables, leads and tubes. If the curve shows a wave form - particularly for the Ne tube - screen the setup better against mains interference, e.g. by grounding the demonstration panel frame.


Fig. 9.2  
Kurve des Hg-Rohres  
Curve of the Hg tube


Fig. 9.3  
Kurve des Ne-Rohres  
Curve of the Ne tube

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The curve form shown on the oscilloscope does not change. Slowly increase the driving potential  $U_1$  and observe the curve: the typical Franck-Hertz curve appears on the screen.

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If the form of the curve is still unsatisfactory even after setting the operating parameters as described in section 3.3.1 / 3, you can optimize it both by tube-specific measures and by matching the parameters to the FH tube you are currently using.

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Increasing the driving potential  $U_1$  dissipates the electron cloud in front of the cathode. As a result, more electrons can be emitted, i.e. the transverse current through the tube is increased.

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When the negative field voltage  $U_3$  is increased, only electrons with a correspondingly high energy can reach the collector, i.e. the collector current is reduced.

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- Maxima are extremely flat, collector current very low:  
increase the driving potential  $U_1$  (or for the Hg tube reduce the oven temperature, see below; on the Ne tube increase the cathode heating voltage, see below).

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- Curve form is not readily apparent:  
Vary the negative field voltage  $U_3$ .

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- Maxima are too high (top of curve is "cut off"): decrease the driving potential  $U_1$  and/or increase the negative field voltage (on the Ne tube, lower the cathode heating voltage, see below).
- Minima approach 0 (bottom of curve is "cut off"): decrease the negative field voltage  $U_3$  and/or increase the driving potential  $U_1$ .

#### 3.3.4.2 Hg tube:

Raising the oven or tube temperature increases the Hg vapor pressure, and thus the probability of electron collisions, i.e. the transverse and collector currents are reduced.

- Collector current abruptly jumps (gas discharge as a result of insufficient oven temperature): immediately set the acceleration voltage  $U_2$  to 0 V and the mode switch (3.1) to RESET. The gas discharge may be observed in the tube as a blue luminance. If the thermal equilibrium was already established when the tube was switched on (if necessary, wait a few minutes), you must increase the oven temperature by changing the set value on the screwdriver potentiometer (2.2).
- Curve is very flat, particularly the first maxima are poorly defined (oven temperature too high): Reduce the oven temperature by changing the set value on the screwdriver potentiometer (2.2). Wait a few minutes until the new thermal equilibrium is reached.

#### 3.3.4.3 Ne tube:

When a fixed acceleration voltage  $U_2$  (approx. 70 V) is supplied, three visible luminance layers should form between the grids when the driving potential  $U_1$  is increased.

Raising the heating voltage increases the temperature of the cathode, and thus the emission current.

- Collector current too high, particularly in oscilloscope observations it is almost impossible to adjust the curve; a bright luminance can be seen between the grids at a low driving potential:  
Lower the cathode temperature by reducing the applied voltage via the screwdriver potentiometer ⑦.
- Collector current too low, maxima barely noticeable; no luminance layers can be seen between the grids even at a high driving potential:  
Raise the cathode temperature by increasing the applied voltage via the screwdriver potentiometer ⑦.

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#### 4 Replacing the fuse

If neither the mains switch ⑥ nor the display (2.6) and the indicator LEDs (3.4) light up, the fuse is blown.

- Pry the insert ① with socket for fuse ② and spare fuse ③ out of its fitting (Fig. 10).
- Replace the old fuse with a new one. Check to make sure that the new fuse has the proper protection values.
- Put in the reserve fuse and replace the insert in the device.


Fig. 10

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#### Important:

If the mains switch lights up, but neither the display nor the LEDs light up, there is a fault in the device which must be repaired at the factory.